


Mepilex® Border Mepilex® Border Lite Product Application Guide


Finger

2" x 5" (5 x 12.5 cm)




Joint

4" x 4" (10 x 10 cm)




Heel

4" x 4" (10 x 10 cm)




Ear

1.6" x 2" (4 x 5 cm)




Toe

2" x 5" (5 x 12.5 cm)



Knee

4" x 12" (10 x 30 cm)



BORDERED FOAM DRESSING WITH SAFETAC® TECHNOLOGY

BENEFITS OF MEPILEX BORDER

- Absorbs moderate to high amounts of exudate
- Self-adherent
- Adheres gently and securely to dry, intact surrounding skin
- Does not adhere to moist wound beds
- Atraumatic to the wound and surrounding skin on removal
- Maintains a moist wound environment
- Minimizes the risk of maceration
- Promotes patient comfort during wear
- May be used under compression bandages
- Shower proof and bacteria proof



BENEFITS OF MEPILEX BORDER LITE

- Conforms well to wound and body contours
- A thin and conformable all-in-one dressing
- Broad selection of sizes
- Minimizes trauma to the wound and pain to the patient during dressing changes
- Shower proof and bacteria proof
- Reduces the risk of maceration
- Does not adhere to moist wound bed
- Adheres gently and securely to dry, intact surrounding skin



WEAR TIME: Up to 7 days

INDICATIONS

Mepilex Border is designed for a wide range of moderate to highly exuding wounds such as:

- Pressure ulcers
- Traumatic wounds
- Secondary healing wounds
- Skin tears
- Leg ulcers
- Foot ulcers
- Partial thickness burns
- Donor Sites

Mepilex Border Ordering Information

Packaged sterile in single packs.

Product Code	Size	Pieces/box	Pieces/case	HCPCS
295200	3" x 3" (7.5 x 7.5 cm)	5	70	A6212
295300	4" x 4" (10 x 10 cm)	5	50	A6212
295400	6" x 6" (15 x 15 cm)	5	50	A6213
295600	6" x 8" (15 x 20 cm)	5	45	A6213

Also Available in  Mepilex® Border Post-Op and  Mepilex® Border Sacrum

WEAR TIME: Up to 7 days

INDICATIONS

Mepilex Border Lite is designed for the management of a wide range of non to low exuding wounds such as:

- Leg and foot ulcers
- Pressure ulcers
- Traumatic wounds, e.g. abrasions, cuts, finger injuries, blisters, and skin tears.
- On compromised and/or fragile skin

Mepilex Border Lite Ordering Information

Packaged sterile in single packs.

Product Code	Size	Pieces/box	Pieces/case	HCPCS
281000	1.6" x 2" (4 x 5 cm)	10	70	—
281100	2" x 5" (5 x 12.5 cm)	5	65	A6212
281200	3" x 3" (7.5 x 7.5 cm)	5	70	A6212
281300	4" x 4" (10 x 10 cm)	5	50	A6212
281500	6" x 6" (15 x 15 cm)	5	50	A6213